

Ergonomic Sonography Chair

Modern, ergonomic design



Description

The ergonomically designed sonography chair provides multiple adjustments for seat height, back height, and backward tilt tension, along with an integrated seat slider allowing for a variety of comfortable postures.

The comfortable medical-grade vinyl with Advanced Beauty Guard provides protection against bacteria.

A foot ring is standard for additional support and positioning.

Specifications

Dimensions:

- **Overall:** 23.5" w x 19.5" d (60 x 50 cm)
- **Seat:**
Width: 97" (49 cm)
Depth: 14.5" - 18.1" (37 cm - 46 cm)
Height: 22.5" - 31.9" (57 cm - 81 cm)
Back: 23.5" w x 17.3" h (60 x 44 cm)
- **Height Adjustable:** 38" - 47.9" (97 cm - 122 cm)
- **Color:** Graphite
- **Weight Capacity:** 250 lb (113.6 kg)
- **Shipping Weight:** 45 lb (20.4 kg)

Ordering Information

058-704 Chair, Ergonomic, Sonography
Includes Foot Ring

mirion.com/biodex