



CERTIFICATE OF COMPLIANCE FOR NUCLEAR MEDICINE SHIPPING SYSTEMS
D.O.T. SPECIFICATIONS 7A, TYPE A

1. Package Identification:

Common Name: Double Dose Pig Shipping System Syringe w/Luer Lock Cap

Catalog Name: PET Shipping System, Double Capped Syringe, Soft Case

Model Number: #001-822 and #001-922 (Retrofit)

2. Preamble:

a) This certificate is issued to satisfy section 178.350 including sections 171.2, 173.403, 173.410, 173.411, 173.412, 173.415, 173.431, 173.465, 173.466, and 178.608 of the Department of Transportation, Hazardous Materials Regulations (49 CFR).

b) The packaging and contents described in item 3 below were found to meet the safety standards set forth in, and tested, according to the more stringent recommendations for general performance testing from:

IATA Dangerous Goods Regulations, 65th Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5.

49 CFR as listed in item 2a above.

c) This certificate does not release the consignor from compliance with any requirement of the regulations of the U.S. Department of Transportation or other applicable regulatory agencies, including the government of any country through or into which the package will be transported.

3. Descriptions of Packaging and Authorized contents, Other Conditions, and References:

a) **Outer Container:** Model #001-922
Sewed Nylon Cloth and Cardboard Shell with metal corner guards.

Dimension: 11.5" x 5.75" x 10.5" with ABS plastic poly foam cushioning and insert support weldment.

b) **Shielded Insert:** Model #001-782-C200
Contour lead shield with open top and bottom and more lead in the vertical center

c) **PET Unit Dose Pig:** Model #001-784 (System accepts 2 pigs)
0.065" polycarbonate with 0.5" lead with single turn twist lock top with 1.44" top lead and 1.62" lead in bottom, and an integral plastic liner.

d) **Absorbent Material:** Model #001-771
3"x3" Quick Solid absorbent material in the liner

e) **Syringe:** B-D plastic 3 ml or 5/6 ml Syringe with Luer



CERTIFICATE OF COMPLIANCE FOR NUCLEAR MEDICINE SHIPPING SYSTEMS
D.O.T. SPECIFICATIONS 7A, TYPE A

Authorized Contents:

Normal Form Radioactive Materials in Type A quantities
Tested with F-18, 235 mCi and 160 mCi in appropriate parts of the shield
as Yellow II. Must comply with A₂ value limits from 49CFR 173.435.

Specifications and Restrictions:

Ground or Air Transportation

4. This certificate is issued on the basis of the test results and safety analysis of the package design.

5. Certification

This document is to certify that to the best of our knowledge, information and belief, the construction methods, the packaging design, and the materials of construction for this specification package comply with the standards described in items 2a, 2b, 2c herein. Any deviation from these conditions must be approved prior to use.

Approved:

Chris Wilczynski
Site Manager, Florham Park • Capintec

Date: 02-14-2024