

Dose Drawing Syringe Shields



Features

- Lightweight, easy-to-use design
- Accommodates standard and safety syringes
- Constructed with a polymer material on each end that acts as a shock absorbers.

Description

Dose Drawing Syringe Shields reduce hand exposure when drawing doses from a shielded vial. Constructed of optically clear high density lead glass, the barrel of the syringe shield offers complete 360° visibility while reducing radiation exposure from Tc-99m by more than 99%.

The shield features quick and smooth syringe insertion with an O-ring seal and anti-roll cap.

Ordering Information

Dose Drawing Syringe Shields:

Accommodates luer and non-luer lock syringes

007-661 Syringe Shield, 3 cc

007-663 Syringe Shield, 5 cc and 6 cc

007-665 Syringe Shield, 10 cc and 12 cc

Dose Drawing Syringe Shields with Replaceable Glass:

Accommodates luer and non-luer lock syringes

007-691 Syringe Shield, 3 cc

007-693 Syringe Shield, 5 cc and 6 cc

007-695 Syringe Shield, 10 cc and 12 cc

Replacement:

127-691 Syringe Shield Replacement Glass, 3 cc

127-693 Syringe Shield Replacement Glass, 5 cc and 6 cc

127-695 Syringe Shield Replacement Glass, 10 cc and 12 cc



Dose Drawing Syringe Shield (glass non-replaceable)

Dose Drawing Syringe Shield with replaceable glass