

Syringe Carrier

Single, Multiple, Small, Large and Shields



Specifications

001-179 Shielded Syringe Carrier, Large

- **Dimensions:**
I.D.: 8.25 x 3 x 2.9 in (21 x 7.6 x 7.4 cm) (l x w x h)
O.D.: 9.5 x 4.4 x 3.5 in (24 x 11.2 x 8.9 cm) (l x w x h)
- **Lead Shielding:**
Sides, top and bottom: .125" thick (.32 cm)
Ends: .25" thick (.64 cm)
- **Weight:** 11.3 lb (5.1 kg)

001-181 Shielded Syringe Carrier, Small

- **Dimensions:**
I.D.: 8 x 1.9 x 1.97 in (20.3 x 4.8 x 5 cm) (l x w x h)
O.D.: 9.25 x 3.4 x 2.6 in (23.5 x 8.6 x 6.6 cm) (l x w x h)
- **Lead Shielding:**
Sides, top and bottom: .125" thick (.32 cm)
Ends: .25" thick (.64 cm)
- **Weight:** 7.5 lb (3.4 kg)

001-182 Shielded Syringe Carrier, Small

- **Dimensions:**
I.D.: 7.5 x 1.7 x 1.7 in (19 x 4.3 x 4.3 cm) (l x w x h)
O.D.: 9.25 x 3.4 x 2.6 in (23.5 x 8.6 x 6.6 cm) (l x w x h)
- **Lead Shielding:**
Sides, top and bottom: .25" thick (.64 cm)
Ends: .5" thick (1.3 cm)
- **Weight:** 11 lb (4.9 kg)

001-180 Shielded Syringe Carrier, Large

- **Dimensions:**
I.D.: 7.7 x 2.6 x 2.7 in (19.6 x 6.6 x 6.9 cm) (l x w x h)
O.D.: 9.5 x 4.4 x 3.5 in (24 x 11.2 x 8.9 cm) (l x w x h)
- **Lead Shielding:**
Sides, top and bottom: .25" thick (.64 cm)
Ends: .5 in thick (1.3 cm)
- **Weight:** 17 lb (7.7 kg)

Description

Shielded Syringe Carriers reduce exposure while storing or transporting radioactive material. The overlapping lid design with snap-latch closure prevents streaming. There are two sizes and thicknesses of lead to choose from. The ends of the carriers are double thick to reduce the exposure from the ends of syringes.

Ordering Information

Shielded Syringe Carrier, .125" lead:

- | | |
|---------|------------------------|
| 001-181 | Syringe Carrier, Small |
| 001-179 | Syringe Carrier, Large |

Shielded Syringe Carrier, .25" lead:

- | | |
|---------|------------------------|
| 001-182 | Syringe Carrier, Small |
| 001-180 | Syringe Carrier, Large |