



Confident. Capable. Connected.

Captus[®] 4000e Thyroid Uptake System

ISO 13485:2016
CERTIFIED



CAPINTEC
A MIRION MEDICAL COMPANY

More Than a Measurement System

Simplify thyroid diagnostics and streamline nuclear medicine workflows with the Captus 4000e. This adaptable system offers modular software, ergonomic design, and seamless EMR integration through DICOM Worklist – backed by the trusted performance of Capintec.



The Captus 4000e is built to support a broad range of nuclear medicine tasks in one compact platform:

- Thyroid Uptake
- Bioassay (I-123, I-125, I-131)
- Wipe Testing with automatic spectrum and trigger alerts
- Automated QA for daily constancy, chi-square, and energy calibration
- Timed Activity Mode and isotope library
- PDF Reporting with full spectrum and graph outputs

Support your I-131 therapies, staff safety, and department compliance with a single versatile system.


Expandable. EMR-Ready.

- DICOM Worklist integration for patient and QC data
- Microsoft Office integration for flexible reporting
- Add-on well counter expands capability to radiation safety
- Shielding options (1 in. or 2 in.) for safety compliance
- Access to parts, service, and support for long-term confidence

Every detail is designed for smoother workflows today – and continued performance tomorrow.



Well Detector (2" Shield)




Neck Phantom

5430-30151	Captus 4000e Thyroid Uptake System
5430-30152	Captus 4000e Thyroid Uptake System with Well
5430-30154	Captus 4000e Thyroid Uptake and Well Counting System with 2" Well Shielding
5430-00007	Captus 4000e DICOM Interface

Related:	
5230-0038	Neck Phantom
5420-2164	Well Detector 2" Shield

mirion.com/captus4000e



Ready to Upgrade Your Workflow?

Recommended for the dependability, ease of usage, and reliability.*

Get expert support, download the spec sheet, or see how the Captus 4000e compares in your department.

Visit mirion.com/captus4000e

*Based on the 2024 Customer Service Satisfaction Survey

Capintec, Inc.

Phone

+1 (800) 631-3826
+1 201 825-9500

Address

7 Vreeland Road,
Florham Park, New Jersey
07932
capintec@mirion.com



CAPINTEC
A MIRION MEDICAL COMPANY

Copyright © 2025 Mirion Technologies, Inc. or its affiliates. All rights reserved. Mirion, the Mirion logo, and other trade names of Mirion products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners.

MKTG-4516 06/2025