

SOUND PRO™ Combination Table

Accommodates ultrasound and echocardiography procedures



Features

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra-wide top with 500 lb patient capacity, suitable for bariatric patients
- Open chassis design and flush-mounted folding side rails for unimpeded patient access
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 15.75" l x 9.5" w
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0 to 80 degrees, designed to accommodate natural body extension and sitting position without slide or shift
- Motorized height, Fowler positioning, Trendelenburg motions and auto-level adjustment
- Naugahyde antimicrobial mattress cover with Advanced Beauty Gard top coat provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller activates table positions
- Allows seated or standing work positions
- Stirrups are designed to cup at the heel for added comfort
- Articulating Scanning Arm Board, adjustable from 0 to 130°
- Foot Controller is ideal for when hands are busy with other functions
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

Optional

- Headrest drops down 15 degrees to allow for carotid and thyroid studies

Description

The Sound Pro Combination Table brings together the imaging features of the Ultra Pro™ and Echo Pro™ tables. Fully equipped for OB/GYN, general ultrasound and echocardiography procedures, the extra-wide tabletop is designed with advanced ergonomics to provide maximum comfort for both patient and sonographer. Designed to conform to natural body extension and sitting position, the motorized Fowler positioning is infinitely adjustable to 80 degrees via hand or foot control. The table is so accommodating, it reduces the total amount of time required to achieve a quality image. Stress patients are off the treadmill and quickly onto the table in imaging position.

The cardiac drop-down cushion releases from either side of the table, providing open access to the patient's left thorax area for an unencumbered apical approach. Flush-mounted side rails enable closer access to the patient, without sonographer contortion. When not in use, rails fold beneath the table, eliminating obstruction.

Patient comfort is ensured, even during lengthy procedures with our Naugahyde antimicrobial mattress cover, which also provides protection against bacteria. Standard accessories available to expand the function of your table include ergonomic stirrups that cup at the heel, an articulating scanning arm board (adjustable from 0 to 130 degrees) and a foot controller for when hands are busy with other functions. The optional headrest is removable and adjustable, providing convenient access and patient comfort during carotid and thyroid procedures.

Biodex | Sound Pro Combination Table



Retractable stirrups expand function and bottom cushion drops out of the way for unobstructed access.

Specifications

• Dimensions:

Overall: 71" l x 30" w (180.3 x 76.2 cm); 35" w (88.9 cm) with side rails

• Weight:

350 lb (159 kg)

• Cardiac Scanning Cutout:

15.75" l x 9.5" w (40 x 24.1 cm) drop-down cushion releases from either side of the table

• Sonographer Cutout:

14.75" l x 8" w (37.4 x 20.3 cm) adjustable sonographer cushion drops down or folds up for maximum access to the patient

• Drop-Down Leg Section:

12.6" l (32 cm)

• Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)

Trendelenburg: 0° to +15°, 0° to -15°

Fowler Positioning: 0° to 80° infinitely adjustable

• Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions

Foot Control: Activates height, Trendelenburg motions and Fowler positioning

• Tabletop:

Three primary sections: torso section features sonographer and patient cutouts with Fowler positioning; center section remains fixed and leg section drops down to 40° and 80° for stirrup access

• Wheels:

5" (12.7 cm) Individual locking swivel casters; central floor-locking system

• Patient Restraints:

One body strap

• Mattress:

Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde antimicrobial mattress cover with Advanced Beauty Gard top coat provides protection against bacteria

• Upholstery Color:

Graphite

• Finish:

Powder coat



Largest cutout in the industry, 15.75" l x 9.5" w.

- **Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating
- **Power:** 115 VAC or 230 VAC
- **Warranty:** Three years parts and labor
- **Certifications:** SGS listed to UL60601-1 and CAN/CSA C22.2 No.601.1.M90 and EN 60601-1 standards
- All Ultrasound Tables are compatible with Hoyer Lifts
- **Primary Function:** General Ultrasound, Echocardiography, Vascular

Ordering Information

058-710 Table, Combination, Sound Pro, 115 VAC
Includes: hand and foot controllers, folding side rails, articulating scanning arm board, retractable stirrups, IV pole and paper dispenser.

058-715 Table, Combination, Sound Pro, 230 VAC
Includes: hand and foot controllers, folding side rails, articulating scanning arm board, retractable stirrups, IV pole and paper dispenser.

Related:

058-704 Chair, Ergonomic, Sonography

Optional:

058-738 Headrest



mirion.com/biodex