

ECONO Echocardiography Table

Extremely versatile,
economically priced



Features

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra-wide top with 500 lb patient capacity, suitable for bariatric patients
- Open chassis design for unimpeded patient access
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 15.75" l x 9.5" w
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0 to 80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height and Fowler positioning
- Naugahyde antimicrobial mattress cover with Advanced Beauty Gard top coat provides protection against bacteria
- Durable hand controller activates table positions
- Allows seated or standing work positions

Optional

- Folding side rails mount flush to table
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest drops down 15 degrees to allow for carotid and thyroid studies

Description

The Econo Echocardiography Table is based on the same platform as the Echo Pro™ Table. It has the same ergonomic features, keeping the safety of the sonographer in mind. A more cost-effective alternative, it is offered without Trendelenburg positioning or the central floor-locking system.

The design includes a cardiac drop-down cushion that provides the largest scanning access of any table, with the ability to release and return from either side of the table. Open access is provided to the patient's left thorax area for an unobstructed apical approach.

The Econo Echocardiography Table has many optional items to further improve sonographer comfort and safety. Choose what best suits your department's specific needs. Redesigned side rails for closer patient access, an articulating scanning arm board, adjustable from 0 to 130 degrees, and a headrest for easier thyroid and carotid access are all available.

Biodex | Econo Echocardiography Table



Adjustable Arm Board 0-130°, locking every 10°.



The cardiac cushion is designed to release from either side of the table.

Specifications

- **Dimensions:**

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails

- **Weight:** 335 lb (152 kg)

- **Cardiac Scanning Cutout:** 15.75" l x 9.5" w (40 x 24.1 cm) drop-down cushion releases from either side of the table

- **Sonographer Cutout:** 14.75" l x 8" w (37.4 x 20.3 cm) adjustable sonographer cushion drops down or folds up for maximum access to the patient

- **Motions:**

Height Adjustable: 23" to 39" (58.4 to 99 cm)

Fowler Positioning: 0° to 80° infinitely adjustable

- **Controls:**

Hand Control: Activates height and Fowler back

- **Tabletop:** Two primary sections: torso section features sonographer and patient cutouts with Fowler positioning; lower body section remains fixed

- **Wheels:** 5" (12.7 cm) Individual locking swivel casters

- **Patient Restraints:** Two body straps

- **Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde antimicrobial mattress cover with Advanced Beauty Gard top coat provides protection against bacteria

- **Upholstery Color:** Graphite

- **Finish:** Powder coat

- **Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating

- **Power:** 115 VAC or 230 VAC

- **Warranty:** Three years parts and labor

- **Certifications:** SGS listed to UL60601-1 and CAN/CSA C22.2 No.601.1.M90 and EN 60601-1 standards

- All Ultrasound Tables are compatible with Hoyer Lifts

- **Primary Function:** Echocardiography

Ordering Information

058-701 Table, Echocardiography, Econo, 115 VAC
Includes hand controller.

058-706 Table, Echocardiography, Econo, 230 VAC
Includes hand controller.

Related:

058-704 Chair, Ergonomic, Sonography

Optional:

058-633 Side Rails, Folding

058-736 Arm Board, Articulating Scanning

058-737 IV Pole

058-611 Dispenser, Paper

058-738 Headrest



mirion.com/biodex