

VASC PRO™ Vascular Ultrasound Table

The ergonomic choice for lower
extremity vascular studies



Features

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra-wide top with 500 lb patient capacity, suitable for bariatric patients
- Reverse Trendelenburg to 30 degrees
- Open chassis design for unimpeded patient access
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0 to 80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height, Fowler positioning, Trendelenburg motions and auto-level adjustment
- Improved foot support with comfortable access to calves without crouching
- Naugahyde antimicrobial mattress cover with Advanced Beauty Gard top coat, provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller
- Allows seated or standing work positions

Optional

- Foot Controller is ideal for when hands are busy with other functions
- Folding side rails mount flush to table
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest drops down 15 degrees to allow for carotid and thyroid studies

Description

The Vasc Pro Vascular Ultrasound Table is designed with ergonomics in mind. It addresses sonographer comfort and safety and proper patient positioning, virtually eliminating repetitive-strain injuries from vascular scanning.

The open chassis design includes an adjustable foot support that accommodates patients of various heights, and optional flush mounted folding side rails to enable close patient access without sonographer contortion.

Fowler positioning to 80 degrees, which accommodates natural body extension and sitting position without slide or shift, an optional Articulating Scanning Arm Board adjustable from 0 to 130 degrees locking every 10 degrees, the thick antimicrobial mattress, and table adjustability all create a comfortable imaging experience for the patient.

Biodex | Vasc Pro Vascular Ultrasound Table



Optional side rails tuck beneath the table when not in use to allow clear access to the tabletop.

Specifications

- **Dimensions:**
 - Overall:** 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails
- **Weight:** 335 lb (152 kg)
- **Motions:**
 - Height Adjustable:** 23" to 39" (58.4 to 99 cm)
 - Trendelenburg:** 0° to +15°, 0° to -30°
 - Fowler Positioning:** 0° to 80° infinitely adjustable
- **Controls:**
 - Hand Control:** Activates height, Trendelenburg motions, Fowler positioning and auto level motions
 - Foot Control (optional):** Activates height, Trendelenburg motions, and Fowler positioning
- **Tabletop:** Two primary sections: torso section features sonographer Fowler positioning; lower body section remains fixed
- **Wheels:** 5" (12.7 cm) Individual locking swivel casters; central floor-locking system
- **Patient Restraints:** Three body straps
- **Foot Support:** Adjustable every 2", up to 14". Stores beneath the table when not in use
- **Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde antimicrobial mattress cover with Advanced Beauty Gard top coat provides protection against bacteria
- **Upholstery Color:** Graphite
- **Finish:** Powder coat
- **Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating
- **Power:** 115 VAC or 230 VAC



The motorized Fowler back has minimal patient shift of no more than 2-3 inches.

- **Warranty:** Three years parts and labor
- **Certifications:** SGS listed to UL60601-1 and CAN/CSA C22.2 No.601.1.M90 and EN 60601-1 standards
- All Ultrasound Tables are compatible with Hoyer Lifts
- **Primary Function:** Vascular

Ordering Information

- | | |
|---------|---------------------------------------------------------------------------------------------|
| 058-732 | Table, Vascular Ultrasound, Vasc Pro, 115 VAC
Includes hand controller and foot support. |
| 058-733 | Table, Vascular Ultrasound, Vasc Pro, 230 VAC
Includes hand controller and foot support. |

Related:

- | | |
|---------|------------------------------|
| 058-704 | Chair, Ergonomic, Sonography |
|---------|------------------------------|

Optional:

- | | |
|---------|----------------------------------|
| 058-741 | Foot Controller |
| 058-633 | Side Rails, Folding |
| 058-736 | Arm Board, Articulating Scanning |
| 058-737 | IV Pole |
| 058-611 | Dispenser, Paper |
| 058-738 | Headrest |



mirion.com/biodex