

CAP-Lift™

Dose Calibrators Accessories

Hands-free dose
measurements



Features

Item Dimensions:

- **Height of lift mechanism:** 19 1/4 in. from top of chamber
- **Diameter of lift stand:** 6 in. at top of chamber
- **Travel:** 8 in. (Standard) Custom travel lengths available for situations where lift stand is not sitting directly on the chamber
- **Control Box:** 4 x 4 x 2 in.



Description

- Reduces hand and body exposure
- Stainless steel construction
- Auto return to elevated dose reception position
- Small size hand control and/or optional foot switch
- Standard and optional custom lengths
- Color coded supply lines
- Perfect alignment with all Capintec™ dose calibrators
- Pneumatic system for faster operations
- No electrical power needed
- Includes foot pedal

Note: This unit requires a low pressure or regulated pressure to operate. A minimum of 25 psi supply line is needed.

Accessories

5120-20302 Replacement Dipper

Ordering Information

5130-30251 Standard CAP-Lift

